

# SeqTrimNext

## Statistics of pre-processing

Plataforma Andaluza de Bioinformática  
Universidad de Málaga

April 23, 2020

# 1 Output Files

SeqTrimNext provides several files, the most interesting ones are in the following directories:

- `output_files`
  - `output.less`, containing an extensive information about the trimming of each sequence. It can be visualised on terminal using the command `less -R`.
  - `used_params.txt`, containing the complete set of parameters used for execution of SeqTrimNext with your data
  - `rejected.txt`, containing a list of rejected sequences together with the reason for their removal.
  - `initial_stats.json`, containing statistics for raw sequences.
  - `stats.json`, containing the statistics of the cleaning process.
  - There is a collection of `folders` that gather sequences with the same MID; each folder contains a `sequences` file (in FASTQ format) with useful reads. There may also exist a file with reads containing low complexity regions. If you want to reconstruct a SFF with the useful segment of each pre-processed read, use `sff_info` file in combination with the original SFF file for the `sfffile` tool.
- `graphs`
  - `size_stats.png`, a graph with the distribution of read lengths in raw data (see Fig. ??).
  - `qualities.png`, a graph to inspect read qualities in raw data (see Fig. ??).
  - `PluginExtractInserts_insert_size.png`, a graph with the distribution of read lengths after SeqTrimNext pre-processing (see Fig. 1).
  - There are other graphs (mostly bar plots) that illustrate the quality of pre-processed reads. All are in PNG format.
- `latex`
  - It is provided as a compressed file `latex.zip` containing all “.tex” files required to compile this document. Graphs are taken from the `graph` folder

# 2 Relevant parameters

In this section, the relevant parameters used in your experiment are shown. Full information about the parameters can be obtained from file `used_params.txt`

## 2.1 General

Plugins applied to every sequence, separated by commas. Order is important

1. PluginIndeterminations
2. PluginFindPolyAt
3. PluginAbAdapters
4. PluginUserContaminants
5. PluginContaminants
6. PluginVectors
7. PluginLowQuality
8. PluginLowComplexity
9. PluginExtractInserts

Remove duplicated (clonal) sequences (using CD-HIT 454)

```
remove_clonality: false
```

Minimum insert size for every trimmed sequence

```
min_insert_size_trimmed: 30
```

Minimum insert size for each end of paired-end reads; true paired-ends have both single-ends longer than this value

```
min_insert_size_paired: 40
```

Seqtrim version

```
seqtrim_version: 2.0.67
```

```
min_sequence_size_raw:
```

## 2.2 Quality

Minimum quality value for every nucleotide

```
min_quality: 20
```

```
window_width:
```

## 2.3 Contaminants

Blast E-value used as cut-off when searching for contaminations

```
blast_evalue_contaminants: 1.0e-10
```

Minimum required identity (%) for a reliable contamination

```
blast_percent_contaminants: 85
```

Minimum hit size (nt) for considering a true contamination

```
min_contam_seq_presence: 40
```

Genus of input data: contaminations belonging to this genus will be ignored

```
genus:
```

Is a contamination considered a source of sequence rejection? (setting to false will only trim contaminated sequences instead of rejecting the complete read)

```

contaminants_reject: true
Path for contaminants database
    contaminants.fasta
    cont_ribosome.fasta

```

### 3 Pre-processing statistics

Next figure is equivalent to Figure ?? but using output reads (useful sequences). The mode is expected to decrease but the shape of the plot should be similar.




Figure 1: Size distribution of the output sequences. Short sequences ( $< \text{min\_insert\_size\_trimmed}$ ) were removed. [PluginExtractInserts\_insert\_size.png]

Summary statistics of the SeqTrimNext analysis. Be careful and read all warnings that are indicating concerns about your data. In the files `initial_stats.json` and `stats.json` can be found a full statistics of your data and SeqTrimNext pre-processing

Input reads:	total	69787796
	Smallest read (bp)	101
	Largest read (bp)	101
	Mode (bp)	0
	Mean (bp)	0.0
Output results:	total	6625528
	Rejected	10590829
	Low complexity reads	107125
	Mode (bp)	91
	Mean (bp)	90.7
	Output paired reads	52464314
	Total output reads	59089842
Linkers:		

Table 1: List of the most frequent Vectors found among your reads

Vectors	sequences
Cloning vector pAAV-MCS, complete sequence.	440624
Cloning vector pKOHPRT complete sequence.	420347
Enterobacteria phage lambda	91825
Cloning vector pVLH/hsp	37902
Cloning vector pWormgate2, complete sequence.	31241

Table 2: List of the most frequent Adapters found among your reads

Adapters	sequences
TruSeq_Universal_Adapter	203694
ABI_Solid3_Adapter_A	189024
ABI_Solid3_GAPDH_Reverse_Primer	141902
Illumina_Single_End_Adapter_1	135816
ABI_Solid3_GAPDH_Forward_Primer	134693

Table 3: List of the most frequent Contaminants found among your reads

Contaminants	sequences
rRNA_long_subunit_Metazoa_Dasytus	78649
Saccharomyces cerevisiae S288c chromosome XII, complete sequence	68497
rRNA_small_subunit_Metazoa_Myotis	62313
Schizosaccharomyces pombe 972h- chromosome III, complete sequence	49903
Aspergillus niger CBS 513.88 clone An03	41155

Table 4: Summary of nucleotides removed in every plugin.

Plugin	Nucleotides	Percent	Warnings
Low Quality	1075001395	Inf %	OK
Low Complexity	52933836	Inf %	OK
Poly T	28807900	Inf %	OK
Poly A	27698228	Inf %	OK
Contaminants	44572920	Inf %	OK
Adapters	58232133	Inf %	OK
Vectors	44920218	Inf %	OK
Indeterminations	174123587	Inf %	OK
Inserts	5370951956	Inf %	iW1

**iW1 Warning!, only Inf % of nucleotides are useful**

## 4 Rejected reads

Input sequences	69787796
Output sequences	6625528
Rejected sequences	10590829
Output paired sequences	52464314
Total output sequences	59089842
Low complexity sequences	107125

Table 5: Summary of reads removed in every plugin.

Case	Number of sequences	Percent	Warnings
	10590829	15.176 %	OK
Short inserts	5933254	8.502 %	rdW2
Empty Inserts	2489841	3.568 %	rdW3
No Valid Inserts	1672614	2.397 %	rdW5
Contaminants	472484	0.677 %	OK
Indeterminations	15490	0.022 %	OK
Low Complexity	7060	0.010 %	OK
Unexpected Vector	86	0.000 %	OK
Total rejected	10590829	15.176 %	OK

**rdW2 Warning!, a 8.502 % of your sequences are too short**

**rdW3 Warning!, a 3.568 % of your sequences are empty (without an insert)**

**rdW5 Warning!, a 2.397 % of your sequences are no valid sequences**

## References

- [1] Falgueras et al. SeqTrim: a high-throughput pipeline for preprocessing any type of sequence reads. *BMC Bioinformatics* 11:38 (2010)
- [2] Weizhong Li & Adam Godzik. Cd-hit: a fast program for clustering and comparing large sets of protein or nucleotide sequences. *Bioinformatics* (2006) 22:1658-9

Thanks you for use SeqTrimNext! Send us any comment to scbi support